



CCC1

**MODEL CCC1 MAXI CASE CART** has two half size shelves for maximum flexibility. Recommended for procedures requiring a larger volume or size of instruments and/or supplies, such as neurosurgery and orthopedic surgery.

Dimensions (body)	Dimensions (overall)	Floor Space	Usable Interior	Shelves	Ship Wt. (lb) est.
46"W	52"W	10.47 sq. ft.	19.43 cu. ft.	2	285
25"D	29"D				
32-7/8"H	40-1/4"H				



CCC2

**MODELS CCC2 AND CCC2E MULTI-PURPOSE CASE CARTS** are the right size for most routine surgical procedures, whether or not your facility utilizes a container system. Model CCC2E has extension shelves for use as a 72" wide back table.

Dimensions (body)	Dimensions (overall)	Floor Space	Usable Interior	Shelves	Ship Wt. (lb) est.
36"W	42"W	8.46 sq. ft.	15.19 cu. ft.	1	235(CCC2)
25"D	29"D		(CCC2/CCC2E)		275(CCC2E)
32-7/8"H	40-1/4"H (CCC2)				
	41-1/4"H (CCC2E)				



CCC3

**MODEL CCC3 MINI CASE CART** is designed for labor and delivery, outpatient surgery or whenever a smaller cart is preferred to transfer instruments and supplies from one area to another.

Dimensions (body)	Dimensions (overall)	Floor Space	Usable Interior	Shelves	Ship Wt. (lb) est.
23-5/8"W	29-5/8"W	5.97 sq. ft.	9.95 cu. ft.	1	175
25"D	29"D				
32-7/8"H	40-1/4"H				





CCC4

**MODEL CCC4 SPACE SAVER CASE CART** is designed to give added capacity by utilizing the cube while minimizing the floor space required by larger carts. Requires the same floor space as the Mini Case Cart CCC3. This cart can handle most routine surgical procedures similar to the Multi-Purpose model CCC2.

Dimensions (body)	Dimensions (overall)	Floor Space	Usable Interior	Shelves	Ship Wt. (lb) est.
23-5/8"W	29-5/8"W	5.97 sq. ft.	14.78 cu. ft.	2	250
25"D	29"D				
47-7/8"H	55-3/8"H				



CCC5

**MODEL CCC5 ULTRA SPACE SAVER CASE CART** offers greater capacity using the same floor space as the Multi-Purpose Cart. It's useful for routine or extensive procedures requiring increased volume/size of instruments and/or supplies.

Dimensions (body)	Dimensions (overall)	Floor Space	Usable Interior	Shelves	Ship Wt. (lb) est.
36"W	42"W	8.46 sq. ft.	22.50 cu. ft.	2	310
25"D	29"D				
47-7/8"H	55-3/8"H				



OCC3

**MODEL OCC3** has an open design that is easy to clean and dry, load and maneuver. Full encircling bumper and solid bottom prevents cross contamination. Ideal when carts go directly to a sterile core, bypassing non-sterile transport.

Dimensions (body)	Dimensions (overall)	Floor Space	Usable Interior	Shelves	Ship Wt. (lb) est.
36"W	42"W	8.46 sq. ft.	N/A	1	200
25"D	29"D				
33-3/16"H	40-1/4"H				



OCC4

**MODEL OCC4** has an open architecture with two wire shelves which allows easy cleaning and drying. Lightweight cart provides economy without sacrificing strength due to solid bottom and encircling bumper.

Dimensions (body)	Dimensions (overall)	Floor Space	Usable Interior	Shelves	Ship Wt. (lb) est.
36"W	42"W	8.46 sq. ft.	N/A	2	175
25"D	29"D				
33-3/16"H	40-1/4"H				



**MODEL 2440 FOLDING CHROME UTILITY CART** is versatile enough for any application requiring transport; i.e. mail, patient records or laboratory materials. Choose any combination of 6" (2441) or 12" (2442) poly-coated wire baskets to complete your cart. The metal frame is designed for easy folding storage. Rubber stoppers support and protect the poly-coated wire baskets. The three-inch, hard rubber casters provide for a quiet ride.

Model No.	Part No.	Features	Dimensions (overall)	Ship Wt.(lb)
2440	2422440000	folding cart	19-1/2"L 29-3/4"W (floor space) 35"H	21
2441	2422441000	6" wire basket	14-1/4"D 24-1/2" W 6"H	5
2442	2422442000	12" wire basket	14-1/2"D 24-3/4"W 12"H	6

**BLICKMAN SUPPLY CARTS** can be used for a variety of applications; as utility carts, general supply carts, par level or exchange carts. Our carts are fabricated from a combination of chrome and stainless. Our carts have a unique design that allows for easy cleanup, protection of supplies and reinforcement of the cart's strength and stability. Its stainless steel bottom pan and wrap around bumper offer maximum protection during transport. Choose the size, style and accessories that meet the needs of your facility. These supply carts come with full-length chrome wire shelves. Seven or fourteen totes can be added to the basic cart for extra storage. Included are a two piece wire mesh back screen, totes with index pockets and label holders to identify supplies and 4 swivel casters with 1 steering lock. Optional accessories are shelf protector, wire shelf dividers, drawer dividers and cart covers.

**GENERAL SUPPLY CARTS (MODELS MEC5167 AND LEC6367)** come either KD or fully assembled with variable shelf adjustment. Designed for easy access, the carts come in two sizes and hold a variety of items on four chrome wire shelves. Wire mesh back screen. Steel chassis with stainless steel bottom and wrap around bumper for maximum protection. Four, six-inch swivel casters, one with swivel lock.

Model No.	Part No.	Dimensions (Overall)	Shelves	Totes	Ship Wt.(lb)
MEC5167	2435267000	51"L 26"W 66"H	4(49"L 21"W)	0	256
LEC6367	2436467000	63"L 26"W 66"H	4(58"L 21"W)	0	282



**EXCHANGE/SUPPLY CARTS (MODELS SETC5467, LETC8067)** come in two sizes, either KD or fully assembled and provide an additional level of organization with seven six-inch deep polyethylene totes in your choice of blue, red or grey. Three chrome wire adjustable shelves and one top shelf, wire mesh back screen and two sizes. Four, six-inch swivel casters, one with swivel lock, allow for easy mobility and smooth steering through confined areas.

Model No.	Part No.	Dimensions (Overall)	Shelves	Totes	Ship Wt.(lb)
SETC5467	2436467100	54"L 26"W 66"H	(3) 39"L 21"W	7	302
LETC8067	2437167000	80"L 26"W 66"H	(3) 46"L 21"W	7	388



For maximum capacity, this **MULTIPURPOSE CART (MODEL METC7467)** utilizes fourteen six-inch deep polyethylene totes, three chrome wire shelves and one top shelf. Wire mesh back screen and four, six-inch swivel casters with one swivel lock for easy mobility and smooth steering.

Model No.	Part No.	Dimensions (Overall)	Shelves	Totes	Ship Wt.(lb)
METC7467	2438367000	74"L 26"W 66"H	(3) 33-1/2"L 21-1/2"W	14	300

## PRODUCT LISTING

### Case Carts

Model No.	Part No.	Dimensions (body)	Dimensions (overall)	Floor Space	Usable Interior	Shelves	Ship Wt. (lb) est.
<b>CCC1</b>	2213331096	46"W 25"D 32-7/8"H	52"W 29"D 40-1/4"H	10.47 sq. ft.	19.43 cu. ft.	2	285
<b>CCC2</b>	2213332096	36"W 25" D 32-7/8"H	42" W 29" D 40-1/4"H	8.46 sq. ft.	15.19 cu. ft.	1	235
<b>CCC2E</b>	2213332596	36"W 25" D 32-7/8"H	42" W 29" D 41-1/4"H	8.46 sq. ft.	15.19 cu. ft.	1	275
<b>CCC3</b>	2213333096	23-5/8"W 25"D 32-7/8"H	29-5/8"W 29"D 40-1/4"H	5.97 sq. ft.	9.95 cu. ft.	1	175
<b>CCC4</b>	2213334096	23-5/8"W 25"D 47-7/8"H	29-5/8"W 29"D 55-3/8"H	5.97 sq. ft.	14.78 cu. ft.	2	250
<b>CCC5</b>	2213335096	36"W 25"D 47-7/8"H	42"W 29"D 55-3/8"H	8.46 sq. ft.	22.50 cu. ft.	2	310
<b>OCC3</b>	2410396000	36"W 25"D 33-3/16"H	42"W 29"D 40-1/4"H	8.46 sq. ft.	N/A	1	200
<b>OCC4</b>	2410496000	36"W 25"D 33-3/16"H	42"W 29"D 40-1/4"H	8.46 sq. ft.	N/A	2	175

## PRODUCT LISTING

### Exchange/Supply and Multi-Purpose Carts

Model No.	Part No.	Description	Dimensions			Ship Wt. (lb) est.
2440	2422440000	Folding Utility Cart	19-1/2"L (49.53)	29-3/4"W (75.57)	35"H (88.90)	21
2441	2422441000	Wire Basket 6"H	14-1/4"D (36.20)	24-1/2"W (62.23)	6"H (15.24)	5*
2442	2422442000	Wire Basket 12"H	14-1/2"D (36.83)	24-3/4"W (62.87)	12"H (30.48)	6*
MEC5167	2435267000	General Supply Cart	51"L (129.54)	26"W (66.04)	66"H (167.64)	256
LEC6367	2436467000	General Supply Cart	63"L (160.02)	26"W (66.04)	66"H (167.64)	282
SETC5467	2436467100	Exchange/Supply Cart	54"L (137.16)	26"W (66.04)	66"H (167.64)	302
LETC8067	2437167000	Exchange/Supply Cart	80"L (203.20)	26"W (66.04)	66"H (167.64)	388
METC7467	2438367000	Multipurpose Supply Cart	74"L (187.96)	26"W (66.04)	66"H (167.64)	300

## OPTIONAL ACCESSORIES FOR EXCHANGE/SUPPLY AND MULTI-PURPOSE CARTS

### Wire Shelf Dividers

Wire Fences	38-1/2" or 45-1/2"
Cart Covers	Various options – consult factory
Divider Boxes (totes)	Various options – consult factory
Dividers	Long or short – for placement inside Divider Boxes
Document Holder (index pocket)	

