

## CHROME



1111

**MODEL 1111 CHROME FOUR-LEG ADJUSTABLE EXAM STOOL** is designed for smooth operation and comfort with a chrome shaft and nylon bushing that rides quietly and cleanly during adjustment. Features a 15"-diameter seat, 3" thick polyfoam padding on durable vinyl upholstery. Stool travels on 2", single-wheel, hard-rubber casters. Choose from 12 upholstery colors.

Model No.	Feature	Part No.	Height Adjustment	Ship Wt. (lb) est.
1111	round foot ring	1021111000	17-1/2"-27-1/2"	15

**MODEL 1113BLK** is our five-leg adjustable exam stool that rests on a black composite base for a contemporary look. Sit in comfort with a 15"-diameter seat, 3" thick polyfoam padding on durable vinyl upholstery. Stool travels on 2", single-wheel, hard-rubber casters. Choose from 12 upholstery colors.

Model No.	Feature	Part No.	Height Adjustment	Ship Wt. (lb) est.
1113BLK	black base	1021113100	20-1/2"-26-1/4"	22
1113WB/BLK	backrest black base	1021113300	20-1/2"-26-1/4"	25



1113BLK



1210WB/BLK

**MODEL 1210WB/BLK FIVE-LEG PNEUMATIC STOOL** offers smooth height adjustment with the touch of the lever on the Supra Pneumatic Lift™. Sit comfortably on the wide 15"-diameter vinyl seat cushion supported by 3" of dense foam padding. The black composite base provides a contemporary look. Two-inch, single-wheel, hard-rubber casters roll easily on carpet or floor. Choose from 12 upholstery colors.

Model No.	Feature	Part No.	Height Adjustment	Ship Wt. (lb) est.
1210BLK	black base backrest	1041210100	22"-29"	25
1210WB/BLK	black base backrest	1041210300	22"-29"	25

## CHROME



**MODELS 1250, 1251 FOOT STOOLS** are constructed of chrome and have 17 weld points to enhance strength and durability. These stools feature a non-slip rubber mat and rubber tips for added comfort. Add the optional handrail – Model 1251.

Model No.	Part No.	Dimensions	Ship Wt. (lb) est.
1250	1011250000	11-3/8"D 14-1/4"W 9"H	8
1251(handrail)	1011251000	11-3/8"D 14-1/4"W 9"H, (35-3/4"HR)	10

*Option Shown: Handrail attachment (model 1251).*



**MODELS 2820 AND 2830 WASTE RECEPTACLES** feature a white painted step-on can with a non-slip pedal and rubber lid guard for secure closure. Black poly pail insert.

Model No.	Feature	Part No.	Dimensions	Ship Wt. (lb) est.
2820	20-Qt.	0832820000	11-1/2"D 11"W 17"H	12
2830	32-Qt.	0832830000	12-1/4"D 11-1/4"W 18-1/4"H	15

## STAINLESS



2840SS

**MODEL 2840SS STAINLESS STEEL WASTE RECEPTACLE** provides a stainless steel step-on can with a quiet, recessed pedal and rubber lid guard for secure closure. Includes a removable poly insert.

Model No.	Feature	Part No.	Dimensions	Ship Wt. (lb) est.
2840SS	20-Qt	0832840000	11-3/4"D 9-3/4"W 17-3/4"H	10

## CHROME

**MODEL 1501 INSTRUMENT STAND** features a five-leg cast aluminum base with wall savers. A Tru- Loc™ friction knob secures the height adjustment of the stainless steel tray. Travels easily on 2", twin-wheel, hooded casters.

Model No.	Feature	Part No.	Height Adjustment	Tray Size	Ship Wt. (lb) est.
1501	Five-leg base	0661501000	29" to 49"	12-5/8"x19-1/8"	12




1501



1510

**DOUBLE-POST AND SINGLE-POST CHROME MAYO STANDS** adjust easily with the Tru- Loc™ friction knob to ensure secure positioning of the stainless steel tray. Travels easily on 2" twin wheel casters (1510) or hard rubber casters (1515).

Model No.	Feature	Part No.	Height Adjustment	Tray Size	Ship Wt. (lb) est.
1510	Double-post	0661510000	34-1/2" to 54"	12-5/8" x 19-1/8"	19
1515	Single-post	0661515000	31-1/2" to 51"	12-5/8" x 19-1/8"	12



## CHROME



2440

**MODELS 2440, 2441, 2442 FOLDING UTILITY CARTS** are versatile enough for any application such as movement of mail, patient records, and laboratory materials. Choose any combination of 6" or 12" deep poly-coated wire baskets to complete your cart. The metal frame is designed for easy assembly and folds for storage. Rubber stoppers support and protect the poly-coated wire racks and provide a quiet ride on four, three-inch hard rubber casters.

Model No.	Feature	Part No.	Dimensions	Ship Wt. (lb) est.
2440	folding cart	2422440000	19-1/2"D 29-3/4"W (floor space) 35" H	21
2441	6" wire basket	2422441000	14-1/4"D 24-1/2"W 6"H	5
2442	12" wire basket	2422442000	14-1/2"D 24-3/4"W 12"H	6

**MODEL 2900 AND 2904 PATIENT PRIVACY SCREENS** are fabricated from sturdy anodized aluminum with stainless steel hinges. The panels are self-extinguishing white vinyl and fold easily for compact storage. All models meet flammability requirements of California Bulletin 117.



2900



2904

Model No.	Feature	Part No.	Dimensions	Ship Wt. (lb) est.
2900	Three-Panel	2822900000	51"L (extended) x 19"W (folded) x 67-1/2"H	9
2904	Three-Panel with casters	2822904000	78-1/4"L (extended) x 26-1/2" W (folded) x 70-1/2" H	19



3400

**MODEL 3400 OVERBED TABLE** moves up or down with a simple one-touch adjustment. Choose an attractive mar-resistant laminate top in wood grain. Security features include a built-in column lock and four-point base. The table travels smoothly on four, two-inch hooded casters.

Model No.	Feature	Part No.	Dimensions	Ship Wt. (lb) est.
3400	wood grain	2833400000	Base-15-1/2"D x 26-1/2"W Top-15"x30"H	27