



7835SS

**Options:** Guard rails and drawers.



7845SS



7821SS

**INSTRUMENT TABLES (HOWARD, SAWYER AND NESTED MODELS)** New and improved design offers complete 18 gauge, 304 stainless construction. Lightweight yet extremely durable. *Sanitary* closed top hem prevents bacterial build-up. *Sealed* reinforcing channels provide structural rigidity. 1-1/4" diameter legs are fully welded for strength and stability. There is a 16" clearance between shelves on the Howard Models. The Sawyer Model has an H-brace in lieu of a second shelf. H-braces are fully welded around the legs with no crevices for dirt to collect. All tables travel on three-inch casters.

Model No.	Part No.	Dimensions (overall)	Feature	Ship Wt. (lb) est.
<b>SHELF MODELS</b>				
7829SS	0117829000	30"L 16"W 34"H	shelf	40
7830SS	0117830000	20"L 16"W 34"H	shelf	32
7831SS	0117831000	33"L 18"W 34"H	shelf	48
7832SS	0117832000	36"L 20"W 34"H	shelf	48
7833SS	0117833000	36"L 24"W 34"H	shelf	60
7834SS	0117834000	48"L 20"W 34"H	shelf	72
7835SS	0117835000	48"L 24"W 34"H	shelf	80
7836SS	0117836000	60"L 24"W 34"H	shelf	112
7837SS	0117837000	72"L 24"W 34"H	shelf	120
3026SS	0113026000	30"W 26"D 34"H	shelf	50
<b>H-BRACE MODELS</b>				
7839SS	0127839000	30"L 16"W 34"H	H-brace	32
7840SS	0127840000	20"L 16"W 34"H	H-brace	24
7841SS	0127841000	33"L 18"W 34"H	H-brace	32
7842SS	0127842000	36"L 20"W 34"H	H-brace	36
7843SS	0127843000	36"L 24"W 34"H	H-brace	40
7844SS	0127844000	48"L 20"W 34"H	H-brace	52
7845SS	0127845000	48"L 24"W 34"H	H-brace	56
7846SS	0127846000	60"L 24"W 34"H	H-brace	96
7847SS	0127847000	72"L 24"W 34"H	H-brace	104
3026SSH	0123026000	30"W 26"D 34"H	H-brace	40

**Options:** Guard rails on 3 sides.

**NESTED INSTRUMENT TABLES** are a space saving alternative with a unique design that requires only 8 sq. ft. of floor space but yields 31 sq. ft. of usable workspace. They can be ordered individually or as a complete set.

Model No.	Part No.	Dimensions (overall)	Ship Wt. (lb) est.
7821SS	0137821000	Complete set of six	300
7822SS	0137822000	48"L 24"W 42"H	110
7823SS	0137823000	44"L 22"W 40"H	90
7824SS	0137824000	40"L 20"W 38"H	80
7825SS	0137825000	36"L 18"W 36"H	60
7826SS	0137826000	32"L 16"W 34"H	50
7827SS	0137827000	28"L 14"W 32"H	40



7850SS



7857SS

**ANESTHESIA/UTILITY TABLES** are popular in the OR and throughout healthcare facilities as a general utility cart. Blickman's superior construction adds a supporting channel underneath the top shelf to enhance its strength; a positive stop to prevent accidental drawer removal and ball bearing slides for an easy-sliding drawer. Travels easily on 3" casters. Choose one to four drawers depending on your requirement.

Model No.	Part No.	Dimensions	Drawers	Ship Wt. (lb) est.
7850SS (Ferguson)	0217850000	20"L 16"W 34"H	1	50
7754SS (George)	0227754000	20"L 16"W 34"H	2	65
7857SS (Saturn)	0237857000	20"L 16"W 34"H	3	70
7859SS (Galaxy)	0247859000	20"L 16"W 34"H	4	110

**Options:** Two 10" drop leaf shelves, crutch tips, larger diameter casters, and guard rail on four sides.



7852SS

**MODEL 7852SS UTILITY TABLE (ALLEN)** is an all purpose table, which offers two solid shelves and a continuous guard rail welded onto the legs on all four sides to protect supplies and equipment. Travels easily on 3" casters.

Model No.	Part No.	Dimensions	Ship Wt. (lb) est.
7852SS	0257852000	20"L 16"W 34"H	30

**Options:** Glides, crutch tips, larger diameter casters.



**MODEL 7856SS UTILITY TABLE (CRESCENT)** is a combination utility-work table, which provides a large work surface and storage in two, four-inch deep drawers housed below the work surface. A continuous guard rail on three sides protects contents. Travels easily on three-inch casters.

Model No.	Part No.	Dimensions	Ship Wt. (lb) est.
7856SS	0317856000	36"L 20"W 34"H	90

**Options:** Glides, crutch tips, larger diameter casters.



**ADJUSTABLE HEIGHT INSTRUMENT TABLES** offer twenty inches of travel to accommodate changes in OR table position during neuro, cardiac and orthopedic surgery. When used in general surgery, the table can be adjusted to accommodate the varying heights of personnel. This model can be adjusted manually from either side with a detachable crank handle. Tables travel easily on four, 3" casters, two casters with brakes.

Model No.	Part No.	Dimensions (overall)	Leg Clearance	Height Adjustment	Ship Wt. (lb) est.
7892SS-0	0157892000	36"L 20"W	28"	36"-56"	140
7893SS-0	0157893000	36"L 24"W	28"	36"-56"	150
7895SS-0	0157895000	48"L 24"W	40"	36"-56"	160
7896SS-0	0157896000	34"L 46"W	26"	36"-56"	175



CSDT6830



CSDT8036 w/ optional casters

**WORK TABLES** are designed for use as prep tables for sterile processing departments but can be used for any application requiring a heavy-duty table with a large surface area. H-brace model comes in three sizes and offers optional shelf placement, 16" high.

Model No.	Part No.	Dimensions (overall)	Feature	Drawers	Ship Wt. (lb) est.
CSDT6830	0336830000	68"L 30"W 36"H	H-Brace	0	240
CSDT6830S1	0336830100	68"L 30"W 36"H	shelf	0	260
CSDT8036	0338036000	80"L 36"W 36"H	H-Brace	0	280
CSDT8036S1	0338036100	80"L 36"W 36"H	shelf	0	300
CSDT9244	0339244000	92"L 44"W 36"H	H-Brace	0	400
CSDT9244S1	0339244100	92"L 44"W 36"H	shelf	0	440

**Options:** Full length overshelf, casters.

Four legs furnished on 68" long tables; all others have six legs.